

Product Datasheet

HSD17B13 Antibody (BF680) (orb1813935)

Description	HSD17B13 Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human
Conjugation	BF680
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human HSD17B13 around the phosphorylation site of Ser33
Target	HSD17B13
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z5P4
Expiration Date	12 months from date of receipt.