

Product Datasheet

PAK6 Antibody (BF680) (orb1813963)

Description	PAK6 Antibody (BF680)
Species/Host	Rabbit
Reactivity	Feline, Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PAK6 around the phosphorylation site of Ser560 RK(p-S)LV
Target	PAK6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NQU5
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.