

Product Datasheet

WHI-P154 (orb181398)

Description	WHI-P154 is a potent JAK3 inhibitor with IC50 of 1.8 μ M, no activity against JAK1 or JAK2, also inh...
Note	For research use only
Formula	$C_{16}H_{14}BrN_3O_3$
MW	376.2
CAS Number	[211555-04-3]
SMILES	<chem>C(O)1=CC=C(NC2=C3C(=NC=N2)C=C(OC)C(OC)=C3)C=C1Br</chem>
Expiration Date	12 months from date of receipt.