


Product Datasheet

DBeQ (orb181411)

Description	DBeQ is a selective, potent, reversible, and ATP-competitive p97 inhibitor with IC50 of 1.5 μM.
Note	For research use only
Formula	C ₂₂ H ₂₀ N ₄
MW	340.421
CAS Number	[177355-84-9]
SMILES	<chem>N([H])(C([H])([H])C1C([H])=C([H])C([H])=C([H])C=1[H])C1C2=C([H])C([H])=C([H])C([H])=C2N=C(N=1)N([H])C([H])([H])C1C([H])=C([H])C([H])=C([H])C=1[H]</chem>
Expiration Date	12 months from date of receipt.



Chemical structure of DBeQ