

## Product Datasheet

### LDN-57444 (orb181419)

<b>Description</b>	LDN-57444 is a reversible, competitive proteasome inhibitor for Uch-L1 with IC50 of 0.88 $\mu$ M, 28-fol...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{11}Cl_3N_2O_3$
<b>MW</b>	397.64
<b>CAS Number</b>	[668467-91-2]
<b>SMILES</b>	<chem>CC(=O)O/N=C\1/C2=C(C=CC(=C2)Cl)N(C1=O)CC3=C(C=CC(=C3)Cl)Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.