

Product Datasheet

MAPT Antibody (AP) (orb1814774)

Description	MAPT Antibody (AP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Tau
Target	MAPT
Preservatives	PBS, Glycerol, rAlbumin
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.