

Product Datasheet

MAPT Antibody (AP) (orb1815499)

Description	MAPT Antibody (AP)
Species/Host	Mouse
Reactivity	Human
Conjugation	AP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human phospho-Tau (Thr231)
Target	MAPT
Preservatives	0.01M PBS(pH7.4)
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.