

Product Datasheet

SARS-CoV-2 Antibody (AP) (orb1815540)

Description	SARS-CoV-2 Antibody (AP)
Species/Host	Mouse
Reactivity	Virus
Conjugation	AP
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Spike S1 Protein (14-685/1213aa)
Target	SARS-CoV-2
Preservatives	0.01M PBS(pH7.4)
Concentration	Lot Dependent
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.