

Product Datasheet

CPV Antibody (AP) (orb1815616)

| | |
|----------------------------|--|
| Description | CPV Antibody (AP) |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | AP |
| Tested Applications | WB |
| Immunogen | Recombinant Canine parvovirus Capsid protein VP2 Fragment C-terminal |
| Target | CPV |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| Uniprot ID | P61826 |
| Expiration Date | 12 months from date of receipt. |