

Product Datasheet

CDV Fusion glycoprotein F0/HRP Antibody (orb1815620)

Description	CDV Fusion glycoprotein F0/HRP Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant CDV Fusion glycoprotein F0 Fragment C-terminal
Target	F
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P12569
Expiration Date	12 months from date of receipt.