

Product Datasheet

ITGAL Antibody (Biotin) (orb1815806)

Description	ITGAL Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ICC
Target	ITGAL
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2b, k
Clonality	Monoclonal
Uniprot ID	P20701
Expiration Date	12 months from date of receipt.