

Product Datasheet

H2-K1 Antibody (PE) (orb1816053)

Description	H2-K1 Antibody (PE)
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	H2-K1
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P04223
Expiration Date	12 months from date of receipt.