

Product Datasheet

ITGA2B Antibody (Biotin) (orb1816086)

Description	ITGA2B Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Target	ITGA2B
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Clonality	Monoclonal
MW	125 kDa
Uniprot ID	P08514
Expiration Date	12 months from date of receipt.