

Product Datasheet

S1 tag Antibody (HRP) (orb1816179)

Description	S1 tag Antibody (HRP)
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide NANNPDWDF
Target	S1 tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.