

Product Datasheet

CPV Antibody (orb1816978)

Description	CPV Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Canine parvovirus Capsid protein VP2 Fragment C-terminal
Target	CPV
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	65 kDa
Uniprot ID	P61826
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.