

Product Datasheet

bs-0294D-BF647 Antibody (BF750) (orb1817042)

Description	bs-0294D-BF647 Antibody (BF750)
Species/Host	Donkey
Reactivity	Goat
Conjugation	BF750
Tested Applications	IF
Immunogen	Native Goat IgG
Target	bs-0294D-BF647
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.