

Product Datasheet

AffiniPure F(ab')₂ Fragment Goat Anti-Rabbit IgG (H+L) (orb1817060)

Description	AffiniPure F(ab') ₂ Fragment Goat Anti-Rabbit IgG (H+L)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Immunogen	Native Rabbit IgG
Preservatives	0.01M PBS (pH7.4)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.