

Product Datasheet

Recombinant human CD71 protein, His (HEK293) (orb1817163)

Description	Recombinant human CD71 protein, His (HEK293)
Target	TFRC
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	89-760/760
Purity	> 90% as determined by SDS-PAGE
MW	76.2 kDa
Uniprot ID	Q62351
Expiration Date	6 months from date of receipt.