

Product Datasheet

Recombinant human LAMP1 protein, His (HEK293) (orb1817167)

Description	Recombinant human LAMP1 protein, His (HEK293)
Endotoxins	Not analyzed
Target	LAMP1
Form/Appearance	Lyophilized or liquid
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	29-382/417
Purity	> 90% as determined by SDS-PAGE
MW	39.7 kDa
Uniprot ID	P11279
Expiration Date	6 months from date of receipt.