

## Product Datasheet

# Recombinant human DR4 protein, C-His-Avi (HEK293) (orb1817171)

<b>Description</b>	Recombinant human DR4 protein, C-His-Avi (HEK293)
<b>Target</b>	TNFRSF10A
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His-Avi
<b>Note</b>	For research use only
<b>Protein Sequence</b>	34-239/468
<b>Purity</b>	> 95% as determined by SDS-PAGE
<b>MW</b>	24.7 kDa
<b>Uniprot ID</b>	<b>000220</b>
<b>Expiration Date</b>	6 months from date of receipt.