

Product Datasheet

Recombinant human M-CSF protein (Active, HEK293) (orb1817190)

Description	Recombinant human M-CSF protein (Active, HEK293)
Endotoxins	< 0.5 EU/μg as determined by LAL.
Target	CSF1
Form/Appearance	Lyophilized powder
Concentration	> 0.5 mg/ml
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. In general: 12 months from date of receipt, -20 to -80°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 3 months, -20 to -80°C under sterile conditions after reconstitution.
Tag	Tag free
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	33-190/554
Purity	> 95% as determined by SDS-PAGE
MW	35-40 kDa
Uniprot ID	P09603
Expiration Date	6 months from date of receipt.