

Product Datasheet

Recombinant human PRL-R protein, C-His (HEK293) (orb1817205)

Description	Recombinant human PRL-R protein, C-His (HEK293)
Endotoxins	Not analyzed
Target	PRLR
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	25-234/622
Purity	> 90% as determined by SDS-PAGE
MW	36 kDa
Uniprot ID	P16471
Expiration Date	6 months from date of receipt.