

## Product Datasheet

# Recombinant human IL-8 protein, His (HEK293) (orb1817209)

<b>Description</b>	Recombinant human IL-8 protein, His (HEK293)
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	CXCL8
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Concentration</b>	> 1mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concent
<b>Protein Sequence</b>	28-99/99
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	9.9 kDa
<b>Uniprot ID</b>	<b>P10145</b>
<b>Expiration Date</b>	6 months from date of receipt.