

Product Datasheet

Recombinant human TLR2 protein, His (orb1817213)

Description	Recombinant human TLR2 protein, His
Endotoxins	Not analyzed
Target	TLR2
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	21-588/784
Purity	> 90% as determined by SDS-PAGE
MW	64 kDa
Uniprot ID	O60603
Expiration Date	6 months from date of receipt.