

## Product Datasheet

# Human TNFRSF17(Tumor Necrosis Factor Receptor Superfamily, Member 17) ELISA Kit (orb1817288)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	0.16-10 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids.
<b>Sensitivity</b>	0.071 ng/mL
<b>Expiration Date</b>	Please enquire.

