

## Product Datasheet

### SA(Sialic Acid) ELISA Kit (orb1817348)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	All
<b>Range</b>	78.13-5000 µg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates and other biological fluids.
<b>Sensitivity</b>	26.1 µg/mL
<b>Expiration Date</b>	Please enquire.

