

## Product Datasheet

# Human CHI3L1(Chitinase-3-Like Protein 1) ELISA Kit (orb1817364)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates and other biological fluids.
<b>Sensitivity</b>	18.1 pg/mL
<b>Expiration Date</b>	Please enquire.

