

Product Datasheet

Rat ACSL4(Acyl-CoA Synthetase Long Chain Family Member 4) ELISA Kit (orb1817375)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Sensitivity	0.061 ng/mL
Expiration Date	Please enquire.

