

## Product Datasheet

# Human LHCGR(Lutropin-choriogonadotropic hormone receptor) ELISA Kit (orb1817382)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.16-10 ng/mL

**Note** For research use only

**Sample Types** serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

**Sensitivity** 0.058 ng/mL

**Expiration Date** Please enquire.

