

## Product Datasheet

### EasyStep Human GH(Growth Hormone) ELISA Kit (orb1817403)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	12.5-800 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma
<b>Sensitivity</b>	0.2 pg/mL
<b>Expiration Date</b>	Please enquire.

