

## Product Datasheet

# EasyStep Human CEA(Carcinoembryonic Antigen) ELISA Kit (orb1817407)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	0.16-10 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma
<b>Uniprot ID</b>	<b>P06731</b>
<b>Sensitivity</b>	7.5 pg/mL
<b>Expiration Date</b>	Please enquire.

