

## Product Datasheet

# Anti-FOLH1 / PSMA Reference Antibody (rosopatamab-MMAE) (orb1817575)

<b>Description</b>	Anti-FOLH1 / PSMA Reference Antibody (rosopatamab-MMAE)
<b>Reactivity</b>	Human, Mouse
<b>Endotoxins</b>	< 0.001EU/μg,determined by LAL method.
<b>Conjugation</b>	MMAE
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Form/Appearance</b>	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Protein Sequence</b>	FOLH1 / PSMA
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	145 Da
<b>Uniprot ID</b>	<a href="#">Q04609</a>
<b>Expiration Date</b>	12 months from date of receipt.