

Product Datasheet

Anti-FOLH1 / PSMA Reference Antibody (rosopatamab tetraxetan) (orb1817576)

Description	Anti-FOLH1 / PSMA Reference Antibody (rosopatamab tetraxetan)
Reactivity	Human, Monkey, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	DM1
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	FOLH1 / PSMA
Isotype	IgG1
Clonality	Monoclonal
MW	145 Da
Uniprot ID	Q04609
Expiration Date	12 months from date of receipt.