

Product Datasheet

Anti-B7-H6 / NCR3LG1 Reference Antibody (Dartmouth patent anti-B7-H6) (orb1817697)

Description	Anti-B7-H6 / NCR3LG1 Reference Antibody (Dartmouth patent anti-B7-H6)
Reactivity	Human, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Protein Sequence	B7-H6 / NCR3LG1
Isotype	IgG
Clonality	Monoclonal
MW	150 Da
Uniprot ID	Q68D85
Expiration Date	12 months from date of receipt.