

Product Datasheet

Anti-PDGFR α / CD140a Reference Antibody (tovetumab) (orb1817809)

Description	Anti-PDGFR α / CD140a Reference Antibody (tovetumab)
Reactivity	Human, Monkey, Mouse, Rat
Endotoxins	< 0.001EU/ μ g, determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative. This antibody is purified through a protein A column.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Protein Sequence	PDGFR α / CD140a
Isotype	IgG1
Clonality	Monoclonal
MW	145.39
Uniprot ID	P16234
Expiration Date	12 months from date of receipt.