

Product Datasheet

Anti-GPNMB Reference Antibody (glembatumumab vedotin) (orb1817869)

Description	Anti-GPNMB Reference Antibody (glembatumumab vedotin)
Reactivity	Human, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	MMAE
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	GPNMB
Isotype	IgG2SA
Clonality	Monoclonal
MW	145.56
Uniprot ID	Q14956
Expiration Date	12 months from date of receipt.