

Product Datasheet

Anti-RHD / CD240d Reference Antibody (roledumab) (orb1817985)

| | |
|----------------------------|--|
| Description | Anti-RHD / CD240d Reference Antibody (roledumab) |
| Reactivity | Monkey |
| Endotoxins | < 0.001EU/μg,determined by LAL method. |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FA, FACS, Kinetics |
| Form/Appearance | Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Protein Sequence | RHD / CD240d |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 147.02 |
| Uniprot ID | Q02161 |
| Expiration Date | 12 months from date of receipt. |