

Product Datasheet

Anti-IL-2Ra / CD25 Reference Antibody (Inolimomab) (orb1818013)

| | |
|----------------------------|--|
| Description | Anti-IL-2Ra / CD25 Reference Antibody (Inolimomab) |
| Reactivity | Human, Monkey |
| Endotoxins | < 0.001EU/μg,determined by LAL method. |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FA, FACS, Kinetics |
| Form/Appearance | Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Protein Sequence | IL-2Ra / CD25 |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 145.1 |
| Uniprot ID | P01589 |
| Expiration Date | 12 months from date of receipt. |