

Product Datasheet

Anti-SIRPa / CD172a Reference Antibody (Hospital for Sick Children patent anti-SIRPA) (orb1818228)

Description	Anti-SIRPa / CD172a Reference Antibody (Hospital for Sick Children patent anti-SIRPA)
Reactivity	Human
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Protein Sequence	SIRPa / CD172a
Isotype	IgG4
Clonality	Monoclonal
MW	145 Da
Uniprot ID	P78324
Expiration Date	12 months from date of receipt.