

## Product Datasheet

# Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab) (orb1818245)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Anti-NRP1 / VEGF165R / CD304 Reference Antibody (vesencumab)   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Endotoxins</b>          | < 0.001EU/μg,determined by LAL method.   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, FA, FACS, Kinetics  |
| <b>Form/Appearance</b>     | Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.       |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only  |
| <b>Protein Sequence</b>    | NRP1 / VEGF165R / CD304  |
| <b>Isotype</b>             | IgG1   |
| <b>Clonality</b>           | Monoclonal   |
| <b>MW</b>                  | 145.22   |
| <b>Uniprot ID</b>          | <b>O14786</b>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |