

## Product Datasheet

# Anti-LILRB2 / ILT4 / CD85d Reference Antibody (JTX-8064) (orb1818262)

<b>Description</b>	Anti-LILRB2 / ILT4 / CD85d Reference Antibody (JTX-8064)
<b>Reactivity</b>	Human, Monkey, Mouse
<b>Endotoxins</b>	< 0.001EU/μg,determined by LAL method.
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Form/Appearance</b>	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LILRB2 / ILT4 / CD85d
<b>Isotype</b>	IgG4
<b>Clonality</b>	Monoclonal
<b>MW</b>	145 Da
<b>Uniprot ID</b>	<a href="#">Q8N423</a>
<b>Expiration Date</b>	12 months from date of receipt.