

Product Datasheet

Anti-PRLR / Prolactin Receptor Reference Antibody (rolinsatamab) (orb1818283)

Description	Anti-PRLR / Prolactin Receptor Reference Antibody (rolinsatamab)
Reactivity	Human, Monkey
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Protein Sequence	PRLR / Prolactin Receptor
Isotype	IgG1
Clonality	Monoclonal
MW	145 Da
Uniprot ID	P16471
Expiration Date	12 months from date of receipt.