

Product Datasheet

Anti-CD79b Reference Antibody (polatuzumab vedotin-piiq) (orb1818307)

Description	Anti-CD79b Reference Antibody (polatuzumab vedotin-piiq)
Reactivity	Human, Monkey
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	MMAE
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	CD79b
Isotype	IgG1
Clonality	Monoclonal
MW	150.53
Uniprot ID	P40259
Expiration Date	12 months from date of receipt.