

Product Datasheet

Anti-LIV-1 / SLC39A6 Reference Antibody ((ladiratuzumab vedotin) (orb1818324)

Description	Anti-LIV-1 / SLC39A6 Reference Antibody ((ladiratuzumab vedotin)
Reactivity	Human, Monkey, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	MMAE
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	LIV-1 / SLC39A6
Isotype	IgG1
Clonality	Monoclonal
MW	146.96
Uniprot ID	Q13433
Expiration Date	12 months from date of receipt.