

## Product Datasheet

# Anti-LOXL2 Reference Antibody (simtuzumab) (orb1818355)

|                                |   |
|--------------------------------|---|
| <b>Description</b>             | Anti-LOXL2 Reference Antibody<br>(simtuzumab)   |
| <b>Reactivity</b>              | Human, Monkey   |
| <b>Endotoxins</b>              | < 0.001EU/μg,determined by LAL<br>method.   |
| <b>Conjugation</b>             | Unconjugated  |
| <b>Tested<br/>Applications</b> | ELISA, FA, FACS, Kinetics   |
| <b>Form/Appearance</b>         | Purified monoclonal antibody supplied<br>in PBS, pH 6.0, without<br>preservative.This antibody is purified<br>through a protein A column.       |
| <b>Concentration</b>           | 1mg/ml  |
| <b>Storage</b>                 | Maintain refrigerated at 2-8°C for up to<br>2 weeks. For long term storage store<br>at -20°C in small aliquots to prevent<br>freeze-thaw cycles |
| <b>Note</b>                    | For research use only   |
| <b>Protein Sequence</b>        | LOXL2   |
| <b>Isotype</b>                 | IgG4  |
| <b>Clonality</b>               | Monoclonal  |
| <b>MW</b>                      | 150 Da  |
| <b>Uniprot ID</b>              | <b>Q9Y4K0</b>   |
| <b>Expiration Date</b>         | 12 months from date of receipt.   |