

## Product Datasheet

# Anti- FcRH5 / IRTA2 / CD307e Reference Antibody (cevastamab) (orb1818357)

<b>Description</b>	Anti- FcRH5 / IRTA2 / CD307e Reference Antibody (cevastamab)
<b>Reactivity</b>	Human, Monkey
<b>Endotoxins</b>	< 0.001EU/μg,determined by LAL method.
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Form/Appearance</b>	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Protein Sequence</b>	FcRH5 / IRTA2 / CD307e
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	145 Da
<b>Uniprot ID</b>	<a href="#">Q96RD9</a>
<b>Expiration Date</b>	12 months from date of receipt.