

Product Datasheet

Anti-ERBB2 / HER2 / CD340 Reference Antibody (trastuzumab deruxtecan) (orb1818442)

Description	Anti-ERBB2 / HER2 / CD340 Reference Antibody (trastuzumab deruxtecan)
Reactivity	Human
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	DXD
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	ERBB2 / HER2 / CD340
Isotype	IgG1
Clonality	Monoclonal
MW	153.75
Uniprot ID	P04626
Expiration Date	12 months from date of receipt.