



Product Datasheet

Caldesmon (Smooth) Cald1 Antibody (Monoclonal, hHCD) (orb18187)

Description	Caldesmon (Smooth) Cald1 Antibody (Monoclonal, hHCD)
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human uterus smooth muscle extract.
Form/Appearance	Lyophilized
Concentration	Add 1ml of PBS buffer will yield a concentration of 100 ug/ml
Storage	At -20°C for one year and at 4°C for one month from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100 ug/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	60584 MW
Uniprot ID	Q62736
Sensitivity	> 5000 cells
Expiration Date	12 months from date of receipt.



IHC(P)
analysis of
Human
Mammary
Cancer ...



IHC(P)
analysis of
Human
Intestinal
Canc...