

Product Datasheet

BMS-986224 (orb1818740)

Description	BMS-986224 is a potent and selective apelin receptor (APJ) agonist with Ki of 0.074 nM, exhibiting s...
Note	For research use only
Formula	$C_{24}H_{23}ClN_4O_6$
MW	498.92
SMILES	<chem>O=C1C(C2=NN=C(CC3=NC=C(Cl)C=C3)O2)=C(O)C(C4=C(OC)C=CC=C4OC)=C(COCC)N1</chem>
Expiration Date	12 months from date of receipt.