

Product Datasheet

BE1218 (orb1818774)

Description	BE1218 is a liver X receptor (LXR) inverse agonist with IC50 values of 9 nM and 7 nM against LXR α a...
Note	For research use only
Formula	C ₃₀ H ₃₀ FNO ₄ S ₂
MW	551.691
CAS Number	[2893967-40-1]
SMILES	<chem>CC1=CC(C)=C(C(C)=C1)S(=O)(N(CC2=C(C=CC=C2)F)CC3=CC=C(C=C3)C4=CC(S(C)(=O)=O)=CC=C4)=O</chem>
Expiration Date	12 months from date of receipt.